

EUROPE

LIEGE SCIENCE PARK • 4102 Seraing • BELGIUM • Tel.: +32 4 372 74 00 • Fax: +32 4 372 75 00
Toll-free: + 800 666 00 123 • info@eurogentec.com • www.eurogentec.com

NORTH AMERICA

ANASPEC - 34801 Campus Drive • Fremont, CA 94555 • USA • Tel.: +1-510-791-9560
Toll-free: +1 800-452-5530 • Fax: +1 510-791-9572 • service@anaspec.com
www.anaspec.com



Technical Data Sheet

Table 4:

3-Step cycling protocol for maximal sensitivity

		Fast Cycling on fast cycler	Regular Cycling
	T°C	Time	
Carry over prevention optional**	50 °C**	2 min.	2 min.
Takyon™ activation	95 °C	3 min.	3 min.
40 Cycles			
Denaturation	95 °C	3 sec.	10 sec.
Annealing	60 °C***	15 sec.	20 sec.
Extension	72 °C	15 sec.	20 - 40 sec.

Technical information

Primer design guidelines

- GC content should be between 30 % and 80 % (ideally 40-60 %).
- Avoid runs of identical nucleotides, especially of 3 or more Gs or Cs at the 3' end.
- The T_m should be between 58 °C and 60 °C.

Custom assay design

The commonly used concentrations for primers are 100 nM. Optimal results may require titration of primers or adjustment of the ratio. The purpose of such a process is to determine the minimum amount of primers required obtaining the most sensitive results with your assay.

Primer titration matrix

Titrate according to the Table 5, perform qPCR and select the concentration which gives the lowest C_q value and clear NTCs. By doing this type of titration it is also possible to compensate for differences up to 2 °C in melt temperature of the primers.

Table 5:

Primer titration matrix

Reverse	Forward		
	50 nM	100 nM	300 nM
50 nM	50 / 50	100 / 50	300 / 50
100 nM	50 / 100	100 / 100	300 / 100
300 nM	50 / 300	300 / 300	300 / 100

MgCl₂ adjustment matrix

Standard MgCl₂ concentration is 2.5 mM but optimal MgCl₂ concentration can vary between assays. If necessary adjust the MgCl₂ concentration with the provided 50 mM MgCl₂ tube. Always prefer optimizing the primer and probe concentrations before the MgCl₂ concentration.

Adjust the amount of water if MgCl₂ is added to the reaction.

For further information please contact our Customer Help Desk:

For Europe:

E-mail: info@eurogentec.com

Tel: +32 4 372 76 65 • Toll-free: + 800 666 00 123

For USA:

E-mail: service@anaspec.com

Tel.: +1-510-791-9560 • Toll-free: +1-800-452-5530

FOR RESEARCH USE ONLY NOTICE TO PURCHASER: LIMITED LICENSE

Licensed dsDNA-Binding Dye Research Kits - Purchase of this product includes an immunity from suit under patents specified in the product insert to use only the amount purchased for the purchaser's own internal research. No other patent rights are conveyed expressly, by implication, or by estoppel. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

UDG - Use of UDG employs U.S. Patents 5,035,996, 5,945,313, 5,683,896 and their foreign counterparts licensed to Eurogentec S.A. from Invitrogen Corporation.

Rox Passive Reference - Purchase of this product includes an immunity from suit to use only the amount purchased for the purchaser's own internal research. No other patent rights are conveyed expressly, by implication, or by estoppel. For information about these rights or on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

Takyon™ is a trademark of Eurogentec.

SYBR® Green I nucleic acid gel stain. This product is provided under a technology license from Life Technologies Corporation to Eurogentec SA. The sale of this product is conditioned on the buyer using the product only in research conducted by the buyer, and not to sell or otherwise commercialize this product or its components (1) in manufacturing; (2) to provide a service, information, or data in exchange for payment; (3) for therapeutic, diagnostic or prophylactic purposes; or (4) for resale, not even if resold for use in research. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@lifetech.com.