

## IFT74 Polyclonal Antibody

<b>Catalog No.</b>	A12672	<b>Category</b>	Polyclonal Antibodies
<b>Applications</b>	WB	<b>Observed MW</b>	80kDa
<b>Cross-reactivity</b>	Human, Mouse, Rat	<b>Calculated MW</b>	42kDa/69kDa

### Immunogen Information

<b>Immunogen</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-372 of human IFT74 (NP_001092694.1).
<b>Gene ID</b>	80173
<b>Swiss prot</b>	Q96LB3
<b>Synonyms</b>	IFT74; BBS20; CCDC2; CMG-1; CMG1; intraflagellar transport 74


### Product information

<b>Source</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification method</b>	Affinity purification
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

### Background

### Recommended Dilutions

WB 1:1000 -  
1:3000



Western blot - IFT74 Polyclonal Antibody (A12672)