

C9orf156 Polyclonal Antibody

Catalog No.	A3434	Category	Polyclonal Antibodies
Applications	WB, IHC	Observed MW	49kDa
Cross-reactivity	Human	Calculated MW	48kDa

Immunogen Information

Immunogen	A synthetic peptide of human C9orf156
Gene ID	51531
Swiss prot	Q9BU70
Synonyms	TRMO; C9orf156; HSPC219; NAP1; tRNA methyltransferase O

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at 4°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000

IHC 1:50 -
1:200



Western blot - C9orf156 Polyclonal Antibody (A3434)