

TP53INP1 Polyclonal Antibody

Catalog No.	A8274	Category	Polyclonal Antibodies
Applications	WB	Observed MW	32kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	18kDa/27kDa

Immunogen Information

Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TP53INP1 (NP_001129205.1).
Gene ID	94241
Swiss prot	Q96A56
Synonyms	TP53INP1; SIP; TP53DINP1; TP53INP1A; TP53INP1B ; Teap; p53DINP1; tumor protein p53-inducible nuclear protein 1


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - TP53INP1 Polyclonal Antibody (A8274)