


Part Number:	PCR-0108LP-TCP-C, PCR-0108LP-TCP-W
Description:	0.1ml 8-strip PCR Tubes and 8-strip optically clear caps

## Specifications

Working volume	< 100μl (0.01ml)
Total well capacity	150μl (0.15ml), Maximum recommended 100μL
Wall thickness	0.0115 inches (+/- 0.005)
Wells	Uniform, thin-walled wells for optimum heat transfer during PCR
Dimensions	 <p>Dimensions in inches, Tolerance +/- 0.0005 in</p>

## Composition

Color	PCR-0108LP-TCP-C (Clear tube) PCR-0108LP-TCP-W (White tube) Caps are optically clear
Material	Medical-grade polypropylene

## Storage and Usability

Storage	Protect from direct sunlight and UV light Store dry at room temperature.
Operational temperature	-20°C to +105°C
Best by date	Five (5) years from date of manufacture
Usage	Designed for single use in PCR applications. Intended for research purposes only.

**Endotoxin (Pyrogen) Testing:** The product is tested to be free of Endotoxins/Pyrogens using the LAL method.

**RNase/DNase Testing:** This product is tested and is free of any detectable RNase/DNase contamination.

**Bovine Spongiform Encephalopathy (BSE) and Transmissible Spongiform Encephalopathy (TSE):** None of these products or raw materials are manufactured from or come into contact with animal material.