

UBR7 Polyclonal Antibody

Catalog No.	A15464	Category	Polyclonal Antibodies
Applications	WB	Observed MW	50kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	47kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 76-425 of human UBR7 (NP_786924.2).
Gene ID	55148
Swiss prot	Q8N806
Synonyms	UBR7; C14orf130; putative E3 ubiquitin-protein ligase UBR7

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background


This gene encodes a UBR box-containing protein that belongs to the E3 ubiquitin ligase family. The protein also contains a plant homeodomain (PHD) in the C-terminus. In mammals, the encoded protein recognizes N-degrons, the destabilizing N-terminal residues of short-lived proteins, which results in ubiquitylation, and proteolysis via the proteasome.

Recommended Dilutions

WB 1:200 -
1:2000



Western blot - UBR7 Polyclonal Antibody (A15464)



Western blot - UBR7 Polyclonal Antibody (A15464)