

NAA50 Polyclonal Antibody

Catalog No.	A4996	Category	Polyclonal Antibodies
Applications	WB, IHC, IF	Observed MW	18kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	9kDa/19kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-169 of human NAA50 (NP_079422.1).
Gene ID	80218
Swiss prot	Q9GZZ1
Synonyms	NAA50; MAK3; NAT13; NAT13P; NAT5; NAT5P; SAN; N-alpha-acetyltransferase 50

Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Background

Recommended Dilutions


WB 1:200 -
1:2000

IHC 1:50 -
1:200


IF 1:50 -
1:200




Western blot - NAA50 Polyclonal Antibody (A4996)




Immunohistochemistry - NAA50 Polyclonal Antibody (A4996)



Immunohistochemistry - NAA50 Polyclonal Antibody (A4996)



Immunohistochemistry - NAA50 Polyclonal Antibody (A4996)



Immunofluorescence - NAA50 Polyclonal Antibody (A4996)