

SNRPA1 Polyclonal Antibody

Catalog No.	A12161	Category	Polyclonal Antibodies
Applications	WB	Observed MW	28kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	28kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human SNRPA1 (NP_003081.2).
Gene ID	6627
Swiss prot	P09661
Synonyms	SNRPA1; Lea1; U2 small nuclear ribonucleoprotein A'


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - SNRPA1 Polyclonal Antibody (A12161)