

HRP-conjugated Mouse anti GST-Tag mAb

Catalog No.	AE027	Category	Tag Antibodies
Applications	WB	Observed MW	
Cross-reactivity		Calculated MW	

Immunogen Information

Immunogen	Recombinant protein of GST-Tag
Gene ID	
Swiss prot	
Synonyms	


Product information

Source	Mouse
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:2000 -
1:10000



Western blot - HRP-conjugated Mouse anti GST-Tag mAb (AE027)