

HIP1R Polyclonal Antibody

Catalog No.	A12026	Category	Polyclonal Antibodies
Applications	WB	Observed MW	130kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	70kDa/119kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-600 of human HIP1R (NP_003950.1).
Gene ID	9026
Swiss prot	O75146
Synonyms	HIP1R; HIP12; HIP3; ILWEQ; huntingtin interacting protein 1 related


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - HIP1R Polyclonal Antibody (A12026)