


Anti-TCFL5 Antibody (A98962)

Specifications:

Name:	Anti-TCFL5 Antibody
Description:	Rabbit polyclonal antibody to TCFL5.
Applications:	WB, ELISA
Recommended Dilutions:	WB: 1:500-1:1000, ELISA: 1:20000
Reactivity:	Human, Mouse
Immunogen:	Synthetic peptide derived from human TCFL5 (amino acids 311-360).
Host:	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugate:	Unconjugated
Purification:	Purified from rabbit serum by antigen affinity chromatography using the immunizing peptide.
Molecular Weight:	48kDa
Product Form:	Liquid
Formulation:	Supplied in Phosphate Buffered Saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, with 150mM NaCl, 0.02% Sodium Azide, and 50% Glycerol.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.


Images:



Western blot analysis of lysates from RAW264.7 cells using Anti-TCFL5 Antibody. The right hand lane represents a negative control, where the antibody is blocked by the immunising peptide.

Anti-TCFL5 Antibody (A98962)

Images continued:



Western blot analysis of RAW264.7 cells using Anti-TCFL5 Antibody.