

## PRODUCT DATA SHEET

## PIG SERUM, STERILE, TRACE-HEMOLYZED, HEAT-INACTIVATED

Product Codes: Pack size: 500 mL

## **Description:**

Whole blood is collected from healthy pigs (mixed breed, either sex) at the time of slaughter, refrigerated within 6 hours, processed within 24 hours and frozen at -20 °C. The serum is later thawed, pooled and passed through a 0.22 micron filter, bottled, and returned to -20 °C. The bottled sterile serum is then thawed in cool water and heat inactivated at 56 °C for 30 minutes. The heat-inactivated serum is then stored at -20 °C.

#### **Physical State:**

Frozen liquid

# Testing:

Hemolysis: ≥ 20 ≤ 50 mg/dL Total Protein: 7.0-10.0 mg /dL

pH: 7.0-8.5

Sterility: Passed, no growth at 14 days

#### Packaging, shipping/storage:

Packaging:

Polypropylene bottle

#### **Storage Temperature:**

-20 °C or below

#### **Shipping Conditions:**

Dry ice

#### **Expiration:**

Product quality is guaranteed to meet Pel-Freez Biologicals' specifications for 1 year from the date of receipt by the customer as long as the product is stored in accordance with the indicated storage conditions.

## **Application Notes**:

Reportedly useful in IHC assays, protein purification.

P/N: 74034 Rev 01