


Anti-EGFR Antibody [A20-E] (A8215)

Specifications:

Name:	Anti-EGFR Antibody [A20-E]
Description:	Rabbit monoclonal (A20-E) antibody to EGFR.
Specificity:	This antibody recognizes the epitope between Ser1120 - Asn1135.
Applications:	IHC-P, IHC-Fr
Recommended Dilutions:	IHC: 1:100-1:200.
Reactivity:	Human
Immunogen:	Peptide derived from C-terminal region of human epidermal growth factor receptor.
Host:	Rabbit
Clonality:	Monoclonal
Clone ID:	A20-E
Isotype:	IgG
Conjugate:	Unconjugated
Product Form:	Liquid
Formulation:	Supplied in Tris Buffered Saline, pH 8.0, with 20 mg/ml BSA and 0.05% Sodium Azide.
Storage:	Store at 4°C.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

Images:



Membranous and cytoplasmic EGFR expression in the colorectal adenocarcinoma. Formalin fixed, paraffin embedded human tissue (4µm section) stained with Anti-EGFR Antibody (A8215).