

## PRODUCT DATA SHEET

### GUINEA PIG COMPLEMENT

**Product Codes**

38002-1  
38002-2

**Pack size:**

10 mL  
20 mL

**Description:**

Serum is processed from whole blood collected from healthy guinea pigs of either sex from US source which are free of visible diseases at the time of blood collection.

**Physical State:**

Frozen liquid

**Testing:**

Hemolytic titer: CH50 > 400 units/mL

**Packaging, shipping/storage:**

**Packaging:**

Polypropylene bottle

**Storage Temperature:**

-70 °C or below

**Shipping Conditions:**

Dry ice

**Expiration:**

Product quality is guaranteed to meet Pel-Freez Biologicals' specifications for 1 year from the date of receipt by the customer as long as the product is stored in accordance with the indicated storage conditions.

**Application Notes:** Reportedly useful in cytotoxicity assays and immunoprecipitation techniques.

P/N: 74077 Rev 01