


Anti-TTK Antibody (A98424)

Specifications:

Name:	Anti-TTK Antibody
Description:	Rabbit polyclonal antibody to TTK.
Applications:	WB, IHC, ELISA
Recommended Dilutions:	WB: 1:500-1:1000, IHC: 1:50-1:100, ELISA: 1:10000
Reactivity:	Human, Mouse
Immunogen:	Synthetic peptide derived from human TTK (amino acids 642-691).
Host:	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugate:	Unconjugated
Purification:	Purified from rabbit serum by antigen affinity chromatography using the immunizing peptide.
Molecular Weight:	97kDa
Product Form:	Liquid
Formulation:	Supplied in Phosphate Buffered Saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, with 150mM NaCl, 0.02% Sodium Azide, and 50% Glycerol.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.


Images:



Immunohistochemical analysis of paraffin-embedded human brain tissue using Anti-TTK Antibody. The right hand panel represents a negative control, where the antibody was pre-incubated with the immunising peptide.

Anti-TTK Antibody (A98424)

Images continued:



Western blot analysis of lysates from HUVEC cells, treated with etoposide 25uM 24H using Anti-TTK Antibody. The right hand lane represents a negative control, where the antibody is blocked by the immunising peptide.