

## ZNF582 Polyclonal Antibody

<b>Catalog No.</b>	A14453	<b>Category</b>	Polyclonal Antibodies
<b>Applications</b>	WB	<b>Observed MW</b>	60kDa
<b>Cross-reactivity</b>	Mouse	<b>Calculated MW</b>	60kDa

### Immunogen Information

<b>Immunogen</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 40-180 of human ZNF582 (NP_653291.1).
<b>Gene ID</b>	147948
<b>Swiss prot</b>	Q96NG8
<b>Synonyms</b>	ZNF582; zinc finger protein 582

### Product information


<b>Source</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification method</b>	Affinity purification
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

### Background

The protein encoded by this gene is a zinc finger protein and putative transcription factor that is highly methylated in cervical cancers. Two transcript variants encoding different isoforms have been found for this gene.

### Recommended Dilutions

WB 1:500 -  
1:2000



Western blot - ZNF582 Polyclonal Antibody (A14453)