

PRODUCT DATA SHEET**HORSE SERUM, DONOR HERD, STERILE, NON-HEMOLYZED****Product Codes:**39321-1
39321-2
39321-3**Pack size:**100 mL
500 mL
3 L**Description:**

Blood is aseptically collected from the donor herd horses (maintained and controlled under the observation of qualified professional personnel) using a closed system method in a designated bleeding lab. The blood is refrigerated within 8 hours of collection and processed within 48 hours prior to storage at -20 °C. Sterilization is achieved using a 0.22 micron filter.

Physical State:

Frozen liquid

Testing:

Hemolysis: ≤ 10 mg/dL
Total Protein: 6.0-9.0 g/dL
pH: 7.0-8.5
Sterility: Passed, no growth at 14 days

Packaging, shipping/storage:**Packaging:**

Polypropylene or PETG bottle

Storage Temperature:

-20 °C or below

Shipping Conditions:

Dry ice

Expiration:

Product quality is guaranteed to meet Pel-Freez Biologicals' specifications for 1 year from the date of receipt by the customer as long as the product is stored in accordance with the indicated storage conditions.

Application Notes:

Reportedly useful in IHC assays as a diluent or blocking agent

P/N: 74033 Rev 01