

MAP1S Polyclonal Antibody

Catalog No.	A17453	Category	Polyclonal Antibodies
Applications	WB	Observed MW	Refer to figures
Cross-reactivity	Human, Mouse, Rat	Calculated MW	109kDa/112kDa

Immunogen Information

Immunogen	A synthetic peptide of human MAP1S.
Gene ID	55201
Swiss prot	Q66K74
Synonyms	MAP1S; BPY2IP1; C19orf5; MAP8; VCY2IP-1; VCY2IP1; microtubule-associated protein 1S

Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:1000

