

MAZ Polyclonal Antibody

Catalog No.	A12921	Category	Polyclonal Antibodies
Applications	WB	Observed MW	60kDa
Cross-reactivity	Human	Calculated MW	17kDa/46kDa/48kDa/51kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-420 of human MAZ (NP_001036004.1).
Gene ID	4150
Swiss prot	P56270
Synonyms	MAZ; PUR1; Pur-1; SAF-1; SAF-2; SAF-3; ZF87; ZNF801; Zif87; myc-associated zinc finger protein


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:1000 -
1:2000



Western blot - MAZ Polyclonal Antibody (A12921)