

AGTRAP Polyclonal Antibody

Catalog No.	A12178	Category	Polyclonal Antibodies
Applications	WB	Observed MW	20kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	11kDa/16kDa/17kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 93-152 of human AGTRAP (NP_001035284.1).
Gene ID	57085
Swiss prot	Q6RW13
Synonyms	AGTRAP; ATRAP; angiotensin II receptor associated protein

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

This gene encodes a transmembrane protein localized to the plasma membrane and perinuclear vesicular structures. The gene product interacts with the angiotensin II type I receptor and negatively regulates angiotensin II signaling. Alternative splicing of this gene generates multiple transcript variants encoding different isoforms.

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - AGTRAP Polyclonal Antibody (A12178)