

RFPL2 Polyclonal Antibody

Catalog No.	A14604	Category	Polyclonal Antibodies
Applications	WB	Observed MW	45kDa
Cross-reactivity	Human, Rat	Calculated MW	32kDa/35kDa/42kDa

Immunogen Information

Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human RFPL2 (NP_006596.2).
Gene ID	10739
Swiss prot	O75678
Synonyms	RFPL2; RNF79; ret finger protein-like 2


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - RFPL2 Polyclonal Antibody (A14604)