

PRODUCT DATA SHEET

RABBIT SERUM, NON-STERILE, TRACE-HEMOLYZED

Product Codes:	<u>Pack size:</u>
31119-1	100 mL
31119-3	3 L
31119-5	500 mL

Description:

Whole blood is collected from young healthy rabbits (Oryctolagus cuniculus) Breeds: New Zealand White or California White that are 8 – 12 weeks of age, either sex and fasted at the time of slaughter, processed, bottled, and stored at -10 °C or below. The serum is then thawed, pooled, and bottled.

Physical State:

Frozen liquid

Testing:

Total Protein: 5.0 - 7.0 g/dLHemoglobin: $20 \le 50 \text{ mg/dL}$ pH: 7.0 - 8.5

Packaging, shipping/storage:

Packaging

Polypropylene bottle with lid

Storage Temperature

-10 °C or below

Shipping Conditions

Dry Ice

Product quality is guaranteed to meet Pel-Freez Biologicals' specifications for 1 year from the date of receipt by the customer as long as the product is stored in accordance with the indicated storage conditions. No expiration date has been assigned.

<u>Application Notes:</u> Reportedly useful in IHC assays, protein purification, and as cell culture and growth media.

P/N: 74006 Rev 03