

OASL Polyclonal Antibody

Catalog No.	A2775	Category	Polyclonal Antibodies
Applications	WB, IHC	Observed MW	63kDa
Cross-reactivity	Human	Calculated MW	29kDa/43kDa/59kDa

Immunogen Information

Immunogen	A synthetic peptide of human OASL
Gene ID	8638
Swiss prot	Q15646
Synonyms	OASL; OASLd; TRIP-14; TRIP14; p59 OASL; p59-OA SL; p59OASL; 2'-5'-oligoadenylate synthetase like

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at 4°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:1000

IHC 1:50 -
1:100



Western blot - OASL Polyclonal Antibody (A2775)