

RBM15 Polyclonal Antibody

Catalog No.	A4936	Category	Polyclonal Antibodies
Applications	WB	Observed MW	120kDa
Cross-reactivity	Human, Mouse	Calculated MW	99kDa/104kDa/106kDa/107kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 530-780 of human RBM15 (NP_001188474.1).
Gene ID	64783
Swiss prot	Q96T37
Synonyms	RBM15; OTT; OTT1; SPEN; RNA binding motif protein 15


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - RBM15 Polyclonal Antibody (A4936)