

ZNF585A Polyclonal Antibody

Catalog No.	A14457	Category	Polyclonal Antibodies
Applications	WB	Observed MW	Refer to Figures
Cross-reactivity	Human	Calculated MW	46kDa/87kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human ZNF585A (NP_001275729.1).
Gene ID	199704
Swiss prot	Q6P3V2
Synonyms	ZNF585A; zinc finger protein 585A

Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000

