

## PRODUCT DATA SHEET

### RABBIT BRAIN CEPHALIN, LIQUID EXTRACT

**Product Codes**41052-1  
41052-2**Pack size:**5 mL  
10 mL**Description:**

Rabbit Brain Cephalin may be used as a phospholipid source in phospholipid-dependent coagulation assays. The phospholipids are isolated from Cholesterol-Free Rabbit Brain Acetone Powder via chloroform extraction, dried, re-suspended 0.85% saline solution, and frozen at -20 °C.

**Physical State:**

Frozen clear liquid

**Testing:**

Phospholipid concentration (PE, PC, PS, PI) reported based on <sup>31</sup>P-NMR analysis.

**Packaging, shipping/storage****Packaging**

Glass Vial

**Storage Temperature**

-20 °C

**Shipping Conditions**

This product is shipped in a frozen state, either with dry ice or in a refrigerated cargo carrier.

**Expiration**

Product quality is guaranteed to meet Pel-Freez Biologicals' specifications for 1 year from the date of receipt by the customer as long as the product is stored in accordance with the indicated storage conditions.

P/N: 74003 Rev 03