

YKT6 Polyclonal Antibody

Catalog No.	A17083	Category	Polyclonal Antibodies
Applications	WB	Observed MW	Refer to figures
Cross-reactivity	Rat	Calculated MW	

Immunogen Information

Immunogen	A synthetic peptide of human YKT6.
Gene ID	10652
Swiss prot	O15498
Synonyms	YKT6; YKT6 v-SNARE homolog

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

This gene product is one of the SNARE recognition molecules implicated in vesicular transport between secretory compartments. It is a membrane associated, isoprenylated protein that functions at the endoplasmic reticulum-Golgi transport step. This protein is highly conserved from yeast to human and can functionally complement the loss of the yeast homolog in the yeast secretory pathway.

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - YKT6 Polyclonal Antibody (A17083)