

PWP2 Polyclonal Antibody

Catalog No.	A5929	Category	Polyclonal Antibodies
Applications	WB	Observed MW	110kDa
Cross-reactivity	Human, Mouse	Calculated MW	102kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 740-919 of human PWP2 (NP_005040.2).
Gene ID	5822
Swiss prot	Q15269
Synonyms	PWP2; EHOC-17; PWP2H; UTP1; periodic tryptophan protein 2 homolog


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:200 -
1:2000



Western blot - PWP2 Polyclonal Antibody (A5929)