

PLCXD2 Polyclonal Antibody

Catalog No.	A14164	Category	Polyclonal Antibodies
Applications	IF	Observed MW	
Cross-reactivity	Human	Calculated MW	34kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-305 of human PLCXD2 (NP_001172035.1).
Gene ID	257068
Swiss prot	Q0VAA5
Synonyms	PLCXD2; PI-PLC X domain-containing protein 2


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

IF 1:50 -
1:100



Immunofluorescence - PLCXD2 Polyclonal Antibody (A14164)