

Apoe Polyclonal Antibody

Catalog No.	A16344	Category	Polyclonal Antibodies
Applications	WB, IF	Observed MW	36kDa
Cross-reactivity	Human, Mouse	Calculated MW	36kDa

Immunogen Information

Immunogen	A synthetic peptide of mouse Apoe
Gene ID	11816
Swiss prot	P02649
Synonyms	APOE; AD2; APO-E; ApoE4; LDLQC5; LPG; AI255918 ; Apo-E; apolipoprotein E

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000

IF 1:50 -
1:200



Western blot - Apoe Polyclonal Antibody (A16344)



Immunofluorescence - Apoe Polyclonal Antibody (A16344)