antibodies

Anti-CDC2 (phospho Thr161) Antibody (A94461)

Specifications:

Name:	Anti-CDC2 (phospho Thr161) Antibody
Description:	Rabbit polyclonal antibody to CDC2 (phospho Thr161).
Specificity:	This antibody detects endogenous levels of CDK1 only when phosphorylated at Thr161.
Applications:	WB, IHC, ELISA
Recommended Dilutions:	WB: 1:500-1:1000, IHC: 1:50-1:100, ELISA: 1:10000
Reactivity:	Human, Mouse, Rat
Immunogen:	Synthetic phospho-peptide derived from human CDC2 around the phosphorylation site of Thr161.
Host:	Rabbit
Clonality:	Polyclonal
lsotype:	lgG
Conjugate:	Unconjugated
Purification:	Purified from rabbit serum by antigen affinity chromatography using the immunizing phospho peptide.
Molecular Weight:	34kDa
Product Form:	Liquid
Formulation:	Supplied in Phosphate Buffered Saline (without Mg2+ and Ca2+), pH 7.4, with 150mM NaCl, 0.02% Sodium Azide, and 50% Glycerol.
Storage:	Shipped at 4° C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

antibodies

Anti-CDC2 (phospho Thr161) Antibody (A94461)

Images:



Western blot analysis of lysates from HeLa cells using Anti-CDC2 (phospho Thr161) Antibody.



ELISA for immunogen phosphopeptide (left) and non-phosphopeptide (right) using Anti-CDC2 (phospho Thr161) Antibody.