

AS3MT Polyclonal Antibody

Catalog No.	A16138	Category	Polyclonal Antibodies
Applications	WB	Observed MW	42kDa
Cross-reactivity	Mouse	Calculated MW	31kDa/41kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 216-375 of human AS3MT (NP_065733.2).
Gene ID	57412
Swiss prot	Q9HBK9
Synonyms	AS3MT; CYT19; arsenite methyltransferase

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

AS3MT catalyzes the transfer of a methyl group from S-adenosyl-L-methionine (AdoMet) to trivalent arsenical and may play a role in arsenic metabolism.

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - AS3MT Polyclonal Antibody (A16138)