


## Anti-TDP43 Antibody (A84627)

### Specifications:

Name:	Anti-TDP43 Antibody
Description:	Goat polyclonal antibody to TDP43.
Applications:	ELISA, WB
Reactivity:	Human
Immunogen:	Synthetic peptide corresponding to Human TARDBP (internal region).
Sequence:	C-HISNAEPKHNSNRQ
Host:	Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugate:	Unconjugated
Purification:	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Concentration:	100 µg at 0.5 mg/ml.
Product Form:	Liquid
Formulation:	Supplied in Tris Buffered Saline, pH 7.30, with 0.02% Sodium Azide and 0.5% BSA.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.


### Images:



Anti-TARDBP Antibody (A84627) (0.3µg/ml) staining of Jurkat lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.

## Anti-TDP43 Antibody (A84627)

Images continued:



Anti-TARDBP Antibody (A84627) (1 $\mu$ g/ml) staining of Human Cerebellum (A), Frontal Cortex (B) and Hippocampus (C) lysates (35 $\mu$ g protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.