


Anti-MYF6 Antibody (A84803)

Specifications:

Name:	Anti-MYF6 Antibody
Description:	Goat polyclonal antibody to MYF6.
Applications:	ELISA, WB
Reactivity:	Human
Immunogen:	Synthetic peptide corresponding to Human MYF6 (internal region).
Sequence:	C-SYRPKQENLEGAD
Host:	Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugate:	Unconjugated
Purification:	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Concentration:	100 µg at 0.5 mg/ml.
Product Form:	Liquid
Formulation:	Supplied in Tris Buffered Saline, pH 7.30, with 0.02% Sodium Azide and 0.5% BSA.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

Images:



HEK293 lysate (10µg protein in RIPA buffer) overexpressing Human MYF6 with MYC tag probed with Anti-MYF6 Antibody (A84803) (0.5µg/ml) in Lane A and probed with anti-MYC Tag (1/1000) in lane C. Mock-transfected HEK293 probed with Anti-MYF6 Antibody (A84803) (0.5mg/ml) in Lane B. Primary incubations were for 1 hour. Detected by chemiluminescence.