

TOMM6 Polyclonal Antibody

Catalog No.	A5240	Category	Polyclonal Antibodies
Applications	WB	Observed MW	13kDa
Cross-reactivity	Mouse	Calculated MW	8kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-74 of human TOMM6 (NP_001127965.1).
Gene ID	100188893
Swiss prot	Q96B49
Synonyms	TOMM6; OBTP; TOM6; translocase of outer mitochondrial membrane 6


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - TOMM6 Polyclonal Antibody (A5240)