PRODUCT INFORMATION



Apigenin 7-O-Glucuronide (hydrate)

Item No. 33754

| Formal Name: | (2S,3S,4S,5R,6S)-3,4,5-trihydroxy-6- ((5-hydroxy-2-(4-hydroxyphenyl)-4- | | | |
|--------------|--|----------|---------------------|----|
| | oxo-4H-chromen-7-vl)oxy)tetrahvdro- | 011 | ОН | |
| | 2H-pyran-2-carboxylic acid, hydrate | OH : | | |
| MF: | $C_{21}H_{18}O_{11} \bullet XH_2O$ | HO | | |
| FW: | 446.4 | | | |
| Purity: | ≥95% HO_ | | | |
| UV/Vis.: | λ _{max} : 270, 335 nm | Π΄ `0΄ U | Ψ U | |
| Supplied as: | A solid | 0 | • XH ₂ O | |
| Storage: | -20°C | | 2 - | ОН |
| Stability: | ≥2 years | | | |
| Item Origin: | Plant/Ixeris sonchifolia | | | |
| | | | 1 110 1 0 | |

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

Apigenin 7-O-glucuronide (hydrate) is supplied as a solid. A stock solution may be made by dissolving the apigenin 7-O-glucuronide (hydrate) in the solvent of choice, which should be purged with an inert gas. Apigenin 7-O-glucuronide (hydrate) is soluble in organic solvents such as DMSO and dimethyl formamide. The solubility of apigenin 7-O-glucuronide (hydrate) in these solvents is approximately 10 mg/ml.

Further dilutions of the stock solution into aqueous buffers or isotonic saline should be made prior to performing biological experiments. Ensure that the residual amount of organic solvent is insignificant, since organic solvents may have physiological effects at low concentrations. Organic solvent-free aqueous solutions of apigenin 7-O-glucuronide (hydrate) can be prepared by directly dissolving the solid in aqueous buffers. The solubility of apigenin 7-O-glucuronide (hydrate) in PBS (pH 7.2) is approximately 1 mg/ml. We do not recommend storing the aqueous solution for more than one day.

Description

Apigenin 7-O-glucuronide is a flavonoid that has been found in J. sigillata and has diverse biological activities¹⁻⁴. It inhibits matrix metalloproteinase-3 (MMP-3), MMP-8, MMP-9, and MMP-13 (IC₅₀s = 12.87, 22.39, 17.52, and 0.27 μ M, respectively).¹⁻⁴ Apigenin 7-O-glucuronide also inhibits protein tyrosine phosphatase 1B (PTP1B), acetylcholinesterase (AChE), and aldose reductase with IC_{50} values of 7.14, 62.96, and 107.1 μ M, respectively.^{2,3} It inhibits LPS-induced production of nitric oxide, prostaglandin E_2 (PGE₂; Item No. 14010), and TNF- α in RAW 264.7 cells when used at a concentration of 100 μ M.⁴ Apigenin 7-O-glucuronide (5 or 10 mg/kg) improves survival in a mouse model of septic shock induced by LPS.

References

- 1. Crasci, L., Basile, L., Panico, A., et al. Planta Med. 83(11), 901-911 (2017).
- 2. Nguyen, D.H., Seo, U.M., Zhao, B.T., et al. Bioorg. Chem. 72, 293-300 (2017).
- 3. Yoo, N.H., Jang, D.S., Yoo, J.L., et al. J. Nat. Prod. 71(4), 713-715 (2008).
- 4. Hu, W., Wang, X., Wu, L., et al. Food Funct. 7(2), 1002-1013 (2016).

WARNING THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFFTY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 04/12/2021

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897 [734] 971-3335 FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.CAYMANCHEM.COM