

ALKBH5 Polyclonal Antibody

Catalog No.	A11684	Category	Polyclonal Antibodies
Applications	WB	Observed MW	51kDa
Cross-reactivity	Human	Calculated MW	43kDa/44kDa/52kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 245-394 of human ALKBH5 (NP_060228.3).
Gene ID	54890
Swiss prot	Q6P6C2
Synonyms	ALKBH5; ABH5; OFOX D; OFOX D1; RNA demethylase A LKBH5


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - ALKBH5 Polyclonal Antibody (A11684)