

[KO Validated] RALBP1 Polyclonal Antibody

Catalog No.	A1140	Category	Polyclonal Antibodies
Applications	WB, IHC	Observed MW	100kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	76kDa


Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 396-655 of human RALBP1 (NP_006779.1).
Gene ID	10928
Swiss prot	Q15311
Synonyms	RALBP1; RIP1; RLIP1; RLIP76; raIA-binding protein 1


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background



Immunohistochemistry - [KO Validated]
RALBP1 Polyclonal Antibody (A1140)




Immunohistochemistry - [KO Validated]
RALBP1 Polyclonal Antibody (A1140)


Recommended Dilutions

WB 1:500 -
1:2000

IHC 1:50 -
1:200



Western blot - [KO Validated] RALBP1
Polyclonal Antibody (A1140)



Western blot - [KO Validated] RALBP1
Polyclonal Antibody (A1140)