

PARP12 Polyclonal Antibody

Catalog No.	A14317	Category	Polyclonal Antibodies
Applications	WB	Observed MW	79kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	79kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 462-701 of human PARP12 (NP_073587.1).
Gene ID	64761
Swiss prot	Q9H0J9
Synonyms	PARP12; ARTD12; MST109; MSTP109; ZC3H1; ZC3HDC 1; poly [ADP-ribose] polymerase 12


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - PARP12 Polyclonal Antibody (A14317)