

In Vivo Buffer pH8.0

In Vivo Buffer pH8.0

Product Details:

ichorbio's *in vivo* buffer pH8.0 is formulated and tested for *in vivo* applications. It is low in endotoxin and has been tested in animal models for toxicity. It can be used with all ichorbio *in vivo* antibodies and biosimilars.

Size:

ichorbio's buffer is available in one standard size: 100ml Please contact us for custom pricing.

Concentration:

1x

Formulation:

50mM Tris, 150mM NaCl. PH 8.0±0.2. Contains no stabilizers or preservatives..

Endotoxin:

? 1.0 EU/mg as determined by the LAL method

Storage:

Store at 2-8°C for up to 12 months.

Use: For research use only (RUO): it is not for diagnostic or therapeutic procedures. **Other buffers:** [In Vivo Buffer 5.0](#)