

ALKBH7 Polyclonal Antibody

Catalog No.	A2331	Category	Polyclonal Antibodies
Applications	WB	Observed MW	26-30kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	24kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-221 of human ALKBH7 (NP_115682.1).
Gene ID	84266
Swiss prot	Q9BT30
Synonyms	ALKBH7; ABH7; SPATA11; UNQ6002; alkB homolog 7


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:200 -
1:2000



Western blot - ALKBH7 Polyclonal Antibody (A2331)