

EZH1 Polyclonal Antibody

Catalog No.	A5818	Category	Polyclonal Antibodies
Applications	WB, IF	Observed MW	85kDa
Cross-reactivity	Human, Mouse	Calculated MW	68kDa/76kDa/80kDa/85kDa/86kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-280 of human EZH1 (NP_001982.2).
Gene ID	2145
Swiss prot	Q92800
Synonyms	EZH1; KMT6B; histone-lysine N-methyltransferase EZH1

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background


Recommended Dilutions

WB 1:500 -
1:2000

IF 1:50 -
1:200



Western blot - EZH1 Polyclonal Antibody (A5818)



Immunofluorescence - EZH1 Polyclonal Antibody (A5818)