

ND5 Polyclonal Antibody

Catalog No.	A12465	Category	Polyclonal Antibodies
Applications	WB, IHC	Observed MW	67kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	68kDa

Immunogen Information

Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of mouse ND5 (NP_904338.1).
Gene ID	17721
Swiss prot	
Synonyms	MT-ND5; MTND5; ND5; NADH dehydrogenase, subunit 5 (complex I)

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background


Recommended Dilutions

WB 1:500 -
1:2000

IHC 1:50 -
1:200



Western blot - ND5 Polyclonal Antibody (A12465)



Immunohistochemistry - ND5 Polyclonal Antibody (A12465)