


Anti-PCSK9 Antibody (A84632)

Specifications:

Name:	Anti-PCSK9 Antibody
Description:	Goat polyclonal antibody to PCSK9.
Applications:	ELISA, WB
Reactivity:	Human
Immunogen:	Synthetic peptide corresponding to Human PCSK9 (internal region).
Sequence:	TRFHRQASKCDSHGT
Host:	Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugate:	Unconjugated
Purification:	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Concentration:	100 µg at 0.5 mg/ml.
Product Form:	Liquid
Formulation:	Supplied in Tris Buffered Saline, pH 7.30, with 0.02% Sodium Azide and 0.5% BSA.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

Images:



HEK293 lysate (10µg protein in RIPA buffer) overexpressing Human PCSK9 with C-terminal MYC tag probed with Anti-PCSK9 Antibody (A84632) (1µg/ml) in Lane A and probed with anti-MYC Tag (1/1000) in lane C. Mock-transfected HEK293 probed with Anti-PCSK9 Antibody (A84632) (1mg/ml) in Lane B. Primary incubations were for 1 hour. Detected by chemiluminescence.