

NTMT1 Polyclonal Antibody

Catalog No.	A7098	Category	Polyclonal Antibodies
Applications	WB, IF	Observed MW	25kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	15kDa/25kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-223 of human NTMT1 (NP_054783.2).
Gene ID	28989
Swiss prot	Q9BV86
Synonyms	NTMT1; AD-003; C9orf32; HONT1A; METTL11A; NRMT ; NRMT1; NTM1A; N-terminal Xaa-Pro-Lys N-methyltransferase 1

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background


Recommended Dilutions

WB 1:500 -
1:2000

IF 1:50 -
1:100



Western blot - NTMT1 Polyclonal Antibody (A7098)



Immunofluorescence - NTMT1 Polyclonal Antibody (A7098)