

STYX Polyclonal Antibody

Catalog No.	A15320	Category	Polyclonal Antibodies
Applications	WB	Observed MW	25kDa
Cross-reactivity	Human	Calculated MW	25kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-223 of human STYX (NP_660294.1).
Gene ID	6815
Swiss prot	Q8WUJ0
Synonyms	STYX; serine/threonine/tyrosine-interacting protein

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

The protein encoded by this gene is a pseudophosphatase, able to bind potential substrates but lacking an active catalytic loop. The encoded protein may be involved in spermiogenesis. Two transcript variants encoding the same protein have been found for these genes.

Recommended Dilutions

WB 1:200 -
1:2000



Western blot - STYX Polyclonal Antibody (A15320)