

ZNF517 Polyclonal Antibody

Catalog No.	A15979	Category	Polyclonal Antibodies
Applications	WB	Observed MW	56kDa
Cross-reactivity	Human	Calculated MW	54kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-180 of human ZNF517 (NP_998770.2).
Gene ID	340385
Swiss prot	Q6ZMY9
Synonyms	ZNF517; zinc finger protein 517


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - ZNF517 Polyclonal Antibody (A15979)