ichorbio

Lusvertikimab Biosimilar - Research Grade

Lusvertikimab Biosimilar - Research Grade

Product Benefits:

ichorbio's lusvertikimab Biosimilar for Research Use Only (RUO)

Target:

IL-7RA

Size:

ichorbio's lusvertikimab biosimilar is available in the following sizes: 1mg, 5mg, 10mg, 20mg, 50mg, 100mg

Isotype:

Human IgG4(S228P/L445P)-Kappa

Host:

CHO Cells

Species Reactivity:

Human

Specificity:

Detects human IL-7RA.

Purification Method:

This monoclonal antibody was purified using multi-step affinity chromatography methods such as Protein A or G depending on the species and isotype.

Concentration:

5mg/ml

Formulation:

Sterile, preservative-free, solution in PBS. BSA and Azide free.

Purity:

>95% by SDS-PAGE

ichor.bio // hello@ichor.bio

Products are for research use only. Not for use in diagnostic or therapeutic procedures. ichorbio, ichorbio logo and all other trademarks are the property of ichorbio LTD © ichorbio 2023 Page 1

ichorbio

Endotoxin:

? 1.0 EU/mg as determined by the LAL method

Storage:

This antibody is stable for at least 4 weeks when stored at 2-8°C. For long term storage, aliquot in working volumes without diluting and store at -20° C or -80° C. Avoid repeated freeze thaw cycles.

Applications:

Functional Assays

Application Notes:

This biosimilar is strictly for research use only (RUO).

ichor.bio // hello@ichor.bio

Products are for research use only. Not for use in diagnostic or therapeutic procedures. ichorbio, ichorbio logo and all other trademarks are the property of ichorbio LTD © ichorbio 2023 Page 2