

## GSTO2 Polyclonal Antibody

<b>Catalog No.</b>	A14156	<b>Category</b>	Polyclonal Antibodies
<b>Applications</b>	IF	<b>Observed MW</b>	
<b>Cross-reactivity</b>	Human	<b>Calculated MW</b>	24kDa/25kDa/28kDa

### Immunogen Information

<b>Immunogen</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 124-243 of human GSTO2 (NP_899062.1).
<b>Gene ID</b>	119391
<b>Swiss prot</b>	Q9H4Y5
<b>Synonyms</b>	GSTO2; GSTO 2-2; bA127L20.1; GSTO2-2; glutathione S-transferase omega-2

### Product information


<b>Source</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification method</b>	Affinity purification
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

### Background

The protein encoded by this gene is an omega class glutathione S-transferase (GST). GSTs are involved in the metabolism of xenobiotics and carcinogens. Four transcript variants encoding different isoforms have been found for this gene.

### Recommended Dilutions

IF 1:50 -  
1:100



Immunofluorescence - GSTO2 Polyclonal Antibody (A14156)