

MYLK4 Polyclonal Antibody

Catalog No.	A10000	Category	Polyclonal Antibodies
Applications	WB	Observed MW	45kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	44kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 119-388 of human MYLK4 (NP_001012418.2).
Gene ID	340156
Swiss prot	Q86YV6
Synonyms	MYLK4; SgK085; myosin light chain kinase family member 4


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:200 -
1:500



Western blot - MYLK4 Polyclonal Antibody (A10000)